BENEfit®

Leading in Orthodontic Implants
psm-medical.com

Quality is never an accident.

It is always the result of high intention, sincere effort, intelligent direction and skillful execution; it represents the wise choice of many alternatives.

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Introduction | Einführung

Skeletal anchorage in orthodontics

An important requirement for successful orthodontic treatment is stable anchorage. Conventional orthodontic anchorage can be limited by periodontal disease and/or missing teeth. Skeletal anchorage has been proven to be very effective. It can be an essential part of non compliance therapy and broadens the options available to the orthodontic practitioner.

Mini Implants are especially advantageous because of the minimally invasive insertion protocol and their relatively low cost.

However, limitations can occur when conventional Mini Implant Systems are used for orthodontic anchorage. In collaboration with leading specialists, **psm** has developed the unique **BENEFit®** TAD System for temporary fixation in the palate. **BENEFit®** offers high precision titanium implant screws and a large variety of pre-fabricated connecting parts such as **BENEplates®**, abutments and palatal expanders.

With the introduction of the **BENEFit®** System **psm** has considerably enhanced the therapeutic diversity of the treatments available for Mesialization, Distalization, RPE as well as uprighting or anchorage of the dentition.

With **BENEFit®**, **psm** was the first company to offer a full range of skeletal anchorage appliances for the palate – a trend that has become „state of the art“ in Orthodontics.

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
Same instruments can be used for several **psm** products. In order to find the right instruments for your use, they are indicated in colour.

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
quattro®

smartlock®

Internal thread for
BENEplate® fixation screw or
the abutment screw.

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Outer double hexagon to prevent
rotation of abutments.

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Optimized thread design leads to easy
insertion without force (cork screw prin-
cipal).

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
BENEFit® Orthodontic Screw 1.5mm, 1ea.
BENEFit® KFO-Schraube 1,5mm, 1Stk.

 All **BENEFit®** screws now also available
 as a pack of 2

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


ST

BENEFit® Orthodontic Screw 2.0mm, 1ea.
BENEFit® KFO-Schraube 2,0mm, 1Stk.

Ø2.3

 Size (Ø x length)
 Abmessung (Ø x Länge)

ST-33-54307	2.3 x 7mm
ST-33-54309	2.3 x 9mm
ST-33-54311	2.3 x 11mm
ST-33-54313	2.3 x 13mm
ST-33-54315	2.3 x 15mm

ST

BENEFit® Orthodontic Screw 2.3mm, 1ea.
BENEFit® KFO-Schraube 2,3mm, 1Stk.
10x

ST-33-54207-10	2.0 x 7mm
ST-33-54208-10	2.0 x 8mm
ST-33-54209-10	2.0 x 9mm
ST-33-54210-10	2.0 x 10mm
ST-33-54211-10	2.0 x 11mm

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


BENEfit® „R“ KFO-Schraube 2,0 mm, 1 Stk.

33-57620 2ea./Stk.

Cylinder for “R” Screw, incl. fixation screws

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Fixierschrauben

BENEplate | BENEplate

33-54400

1ea. / 1Stk.

BENEplate, L, incl. 2 fixation screws, stainless steel

BENEplate, L, inkl. 2 Fixierschrauben, Implantatstahl

33-54402

1ea. / 1Stk.

BENEplate, S, incl. 2 fixation screws, stainless steel

BENEplate, S, inkl. 2 Fixierschrauben, Implantatstahl

33-54424

1ea. / 1Stk.

BENEplate, XS, incl. 2 fixation screws, stainless steel

BENEplate, XS, inkl. 2 Fixierschrauben, Implantatstahl

33-54429

1ea. / 1Stk.

BENEplate, L, median, incl. 2 fixation screws, with 1.1mm wire (12cm), stainless steel

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33-54409

1ea. / 1Stk.

BENEplate, S, median, incl. 2 fixation screws, with 1.1mm wire (12cm), stainless steel

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33-54423

1ea. / 1Stk.

BENEplate, XS, median, incl. 2 fixation screws, with 1.1mm wire (12cm), stainless steel

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33-54428

1ea. / 1Stk.

BENEplate, L, median, incl. 2 fixation screws, with 0.8mm wire (12cm), stainless steel

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
33-54408

1ea. / 1Stk.

BENEplate, S, median, incl. 2 fixation screws, with 0.8mm wire (12cm), stainless steel

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33-54475

1ea. / 1Stk.

BENEplate, S, median, incl. 2 fixation screws, with 0.8 + 1.1mm wire (12cm), stainless steel

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33-54477

1ea. / 1Stk.

BENEplate, S, median, incl. 2 fixation screws, with bilateral 1.1mm wire (6cm), stainless steel

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BENEplate | BENEplate

33-54472 1ea. / 1Stk.

BENEplate, L, paramedian, inkl. 2 fixation screws, with 1.1mm wire (12cm), stainless steel

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33-54471 1ea. / 1Stk.

BENEplate, S, paramedian, inkl. 2 fixation screws, with 1.1mm wire (12cm), stainless steel

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33-54422 1ea. / 1Stk.

BENEplate, XS, paramedian, inkl. 2 fixation screws, with 1.1mm wire (12cm), stainless steel

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33-54474 1ea. / 1Stk.

BENEplate, S, paramedian, inkl. 2 fixation screws, with 0.8mm wire (12cm), stainless steel

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

33-54480

1ea. / 1Stk.

BENEplate, L, median, incl. 2 fixation screws, with central 1.1mm wire (12cm), stainless steel, acc. Dr. Nienkemper, Düsseldorf, Germany

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
33-54485

1ea. / 1Stk.

BENEplate, L, T-Mesialslider plate, median, incl. 2 fixation screws, with central 1.1mm wire + 0.8mm wire anteriorly (12cm), stainless steel


BENEplate, L, T-Mesialslider plate, median, inkl. 2 Fixierschrauben, mit 1,1mm mittigem Draht + 0,8mm Draht anterior (12cm), Implantatstahl

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

A device for molar distalization. The upper molars will be derotated, uprighted and distalized in only two or three months. By converting the system in a transpalatal bar, it can be used as an anchor for the molars. The anterior teeth can be retracted and aligned simply.

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**33-54465**

1ea. / 1Stk.

Hyrax® Ring, incl. 1 fixation screw, with 1.1mm wire (6cm) laterally, stainless steel

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
**33-54461**

1ea. / 1Stk.

Hyrax® Ring, incl. 1 fixation screw, with 1.1mm wire (12cm), stainless steel

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
**33-54464**

1ea. / 1Stk.

Hyrax® Ring, incl. 2 fixation screws, with 1.1mm wire (12cm), stainless steel

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**33-54750**

1ea. / 1Stk.

Hyrax® Ring, incl. 1 fixation screw, with 1.1mm wire (6cm), stainless steel

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BENEslider | BENEslider – Distalization

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
33-54544

2ea. / 2Stk.

BENEtube, with wire 1.1mm, 4cm, acc. to Dr. T. Banach

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

The BENEslider

The most classical indication for skeletal anchorage using Mini-Implants with abutments is molar-distalization and/ mesialization in the maxilla. By means of the so called "**BENEslider**" the molars can be bodily distalized and/or mesialized using the Mesialslider. In many cases, tooth extraction or dental implants can be avoided.

To couple two **BENEfit®** Mini-Implants, a **BENEplate** with prewelded wire is connected to the implants and the **BENEslider** parts (mobilizer, springs, **BENETube**) are placed on the wire. Active force is applied by springs (240g or 500g) and activated with the Mobilizer.

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33-57620

2ea. / 2Stk.

Cylinder for "**R**" **Screw**, incl. fixation screws

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33-54460

1ea. / 1Stk.

Abutment Standard with 1.1mm wire (12cm), stainless steel

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33-54445

1ea. / 1Stk.

Abutment with loop, acc. to Dr. S. Renger, St. Louis, France

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33-54466

1ea. / 1Stk.

BENEFit® PEEK Abutment, incl. 1 fixation screw

BENEFit® PEEK Abutment, inkl. 1 Fixierschraube



33-54469

1ea. / 1Stk.

BENEFit® PEEK Abutment, 15° angled, incl. 1 fixation screw

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Hybrid Hyrax® RPE* | Hybrid Hyrax® GNE*


The **Hyrax®** Ring and **Hyrax®** Abutment facilitate the construction and fixation of a **Hybrid Hyrax®** appliance. It allows welding or soldering with the expansion screw as well as little rotational movements when connecting the appliance to the **BENEFIT®** Screws.

The **Hyrax®** Ring is supplied with 1 fixation screw each.

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33-54465

1ea. / 1Stk.

Hyrax® Ring, incl. 1 fixation screw,
with 1.1mm wire (6cm) laterally, stain-
less steel

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33-54463

1ea. / 1Stk.

BENEFIT® Hyrax®, Abutment
BENEFIT® Hyrax®, Abutment

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bmx® – BENEFIT® Maxillary Xpander | bmx® – BENEFIT® Maxillary Xpander

After using the **BENEFIT®** Screw for the **Hybrid Hyrax®** appliance very successfully, in the last years the **bmx®** now follows the Orthodontists request for an entirely **BENEFIT®** mini-implant anchored appliance for rapid maxillary expansion on young patients. Using the De Franco parallel guide, (Page 35) the distance and the parallelity of the **BENEFIT®** Screws can be exactly determined, on the patient directly. An impression is no longer necessary. Alternatively, visualisation can, of course, be made by means of digital imaging or traditional impression taking.

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bmx®, palatal expander


bmx®, Gaumennahtspreizer

MARPE – Miniscrew Assisted Rapid Palatal Expansion
smartlock®

 [Video](#)

Quadexpander

BENEFIT® DIRECT Screws | BENEFIT® DIRECT Schrauben
WHY Smartlock BENEFIT DIRECT screws for your perfect MARPE result?

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
BENEFIT® DIRECT Screw 2.0 mm
BENEFIT® DIRECT Schraube 2,0 mm

33-64410

2ea. / 2 St.

Hyrax® Ring DIRECT 1-hole, stainless steel

Hyrax® Ring DIRECT 1-Loch, Implantatstahl

2X
ST
ST-33-64209-2 2.0 x 9mm

ST-33-64211-2 2.0 x 11mm

BENEFIT® DIRECT Screws in a package of 2 ea.

BENEFIT® DIRECT Schrauben im VPE 2 St.

33-64414


1ea. / 1 St.

Hyrax® Ring DIRECT 1-hole with wire 1.1mm, 6cm, stainless steel

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Mouse Trap – for Molar Intrusion | Mausefalle – zur Molarenintrusion

consisting of | bestehend aus:

**ST-33-54209-2**

Page | Seite 8

BENEFit® Screw 2.0 x 9mm, 2ea.**BENEFit®** Screw 2,0 x 9mm, 2Stk.**33-54408**

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Scanbodies | Scanbodies

As a further step into digitalization of the **BENEFIT®** System, these scanbodies were developed in cooperation with Orthodontist Simon Graf (Belp, Switzerland) to ensure visualization of a detailed rendering.

The Scanbody Standard is used right after implant insertion. The Scanbody mini can be screwed onto the implants even after several months of treatment.

Developed in cooperation with:
Simon Graf, Orthodontist, Belp, Switzerland

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**33-10760**

Scanbody PEEK Standard
Scanbody PEEK Standard

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10 Stück/Pkg.

Sterile packaging | Sterile Verpackungen

Implants – in single sterile packaging

The outer and inner packaging contains labels indicating UDI codelike article number, LOT number, expiration date, etc. for the patient folder and for full traceability.

The sterile screws are contained in a small colour-coded screw carrier which protects the screw and at the same time acts as a pick up tray. This carrier is placed in a clear, sealed blister package. The screw can easily be removed from the carrier by means of the self retaining screwdriver.

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Note | Hinweis:

The colors of the screw carrier and the correct drills are matching.

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Instruments Instrumente

Drills | Bohrer

The instruments can be used for several **psm** products. In order to find the right instruments for your use, they are indicated in color.

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10-67513

1ea. / 1Stk.

Drill, 1.4x33mm WL 15mm, red, for 2.0mm screws
Bohrer, 1,4x33mm AT 15mm, rot, für 2,0mm Schrauben

11-18452

1ea. / 1Stk.

Drill, 1.8x33mm WL 15mm, grey, for 2.3mm screws
Bohrer, 1,8x33mm AT 15mm, grau, für 2,3mm Schrauben

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The extension of 16 mm works for all rotary instruments of the **quattro®** and **BENEFIT®**-System with a latch type DENTAL attachment (max. rpm 600).

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
Drill stop | Bohrerstop

The drill stop prevents the clinician from unintentionally drilling too deep. They assure an accurate working length from 7 to 13 mm in 2 mm increments. By sliding the stop sleeve over the appropriate **psm** drill, the working length is adjusted to the desired depth. Removal of the sleeve is made easy by using the buffer.

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**33-18307**


Drill-Stop sleeve, 7mm
Bohrer-Stophülse, 7mm

1ea. / 1Stk.

**33-18309**

Drill-Stop sleeve, 9mm
Bohrer-Stophülse, 9mm

1ea. / 1Stk.

**33-18311**

Drill-Stop sleeve, 11mm
Bohrer-Stophülse, 11mm

1ea. / 1Stk.

**33-18313**

Drill-Stop sleeve, 13mm
Bohrer-Stophülse, 13mm

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33-18288

1ea. / 1Stk.

Laboratory Screwdriver, fixed
Labor Schraubendreher, fest

33-54519

1ea. / 1Stk.

Activating Screwdriver for mobilizer, not autoclavable, acc. Dr. Michael Schubert, Regensburg, Germany

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33-10903

1ea. / 1St.

BENEFIT® / Digital Motion Screw holder with stop

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BENEFIT® DIRECT TX Screw holder, 28 mm

BENEFIT® DIRECT TX Schraubenhalter, 28 mm

BENEFIT® DIRECT Screw | BENEFIT® DIRECT Schraube

Battery Screwdriver – iSD900 | Akku Schraubendreher iSD900

Cordless Battery Screwdriver with Torque Calibration System for the controlled fixation of psm Miniscrews.

The iSD 900 will guarantee precise torque control due to the unique – easy to use – torque calibration system (10 – 40 Ncm). Time-saving and precise, riskfree insertion and tightening of screws and abutments. The complete set contains a reliable and user friendly handpiece for Dental attachments. Reverse mode is indicated with sound. Standard rechargeable AAA Ni-Mh batteries can be used.

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
Features | Vorteile

- 10 - 40 Ncm torque in 1 or 5 increments
- 5, 10, 15, 25 rpm turning speed
- Reverse mode with sound
- Calibration facility
- Standard rechargeable batteries

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Mobilholder | Mobilholder

The new Mobilholder is intended to move the **BENEFit®** mobilizers into the final position prior to tightening the set screw for distalization (fig. 2, 3) or mesialization (fig. 1).

developed in cooperation with ZTM W. Bervoets

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
Mobilholder
acc. to ZTM W. Bervoets
Mobilholder
nach ZTM W. Bervoets

Handles | Handgriffe

The manually turned unit is connected to your existing contra-angled handpiece (blue 1:1) that enables its use like an angled screwdriver.

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The soft tissue probes are infinitely adjustable and are fixed by means of the normal **BENEFIT®** instruments (i. e. **33-18288**) with set screws in the desired position.

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33-18453

1ea. / 1Stk.

Connecting rod for soft tissue probe

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33-18452

2ea. / 2Stk.

Soft Tissue Probe

Mucosa Markierungssonde

33-18450

1ea. / 1Stk.

BENEFIT® Screw Location Probe Set complete, 4 parts

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The screw location probe acc. to Dr. Schlichter identifies the target points for the median or paramedian placement of the **BENEFIT®** Screws. The precise positioning of the **BENEFIT®** Implants improves the location and parallelism with each other. It also increases stability and simplifies the chairside bending of the plate. Additionally the lab time is reduced due to less bending.

Developed in cooperation with Dr. Kai Schlichter, Dres. Werling & colleagues, Offenbach a. d. Q., Germany.

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


33-10912 12mm

BENEsleeves double, parallel, 2 ea.
BENEsleeves doppel, parallel, 2 Stk.

The matching drill (**33-10914**) and screwholder with stop (**33-10903**) you will find on page 35.

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The screwholder (33-10903) and drill (33-10914) can be used with digital motion (page 34), De Franco parallel guide and individually made insertion guides (i.e. TADMAN).

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95-13012 1ea. / 1Stk.

Mouse Trap Intrusion
Mausefalle Intrusion

95-13013 1ea. / 1Stk.

Hybrid Hyrax® Distalizer
Hybrid Hyrax® Distalizer

95-13014 1ea. / 1Stk.

T-Mesial/Distalslider
T-Mesial/Distalslider

95-13001-15 1ea. / 1Stk.

BENEslider® mesial/distal
BENEslider® mesial/distal

V-90-994-00 1ea. / 1Stk.

BENEFit® Model Case for
9 models, empty
BENEFit® Modellkoffer für
9 Modelle, leer

V-90-994-02 1ea. / 1Stk.

BENEFit® Model Case for
18 models, empty
BENEFit® Modellkoffer für
18 Modelle, leer

median

paramedian

95-13017-15 1ea. / 1Stk.

BENEslider® mesial/distal
BENEslider® mesial/distal

Color | Farbe

Please select the desired **BENEFit**® model/s and choose the color from the QR code file.

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Catalogues | Kataloge

More innovative **psm** products are presented in our other catalogues. Please request it from your local **psm** distribution partner or download from our website www.psm-medical.com.

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KAT-034

smartlock

The new **psm – smartlock®** system is largely based on the successful and proven **BENEFIT®** system. **smartlock®** focuses on the essentials: Unique comfort as well as safe and easy handling that doctors and patients can rely on. The trend-setting concept, combined with the best materials and the most modern manufacturing processes, guarantees a treatment result at the highest level.

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
KAT-003

quattro

The **quattro®** Mini Anchor System is known as one of the first and most successful orthodontic Mini-Implant Systems worldwide. Based on this system the **BE-NEFIT®** screws were developed. This enables the user to utilize the entire scope of orthodontic implantation with unmatched possibilities and compatible instrumentation.

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KAT-010

mentoplate

The **mentoplate®** system provides a perfect solution for enhanced stability with special orthodontic bone anchored plates for class III correction. The System includes Titanium plates which are fixed to the bony surface with selfdrilling screws.

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EASY DRIVER® THE GUIDING SYSTEM FOR TADS

Digital skeletal anchorage | Digitale skelettale Verankerung

A new and innovative method for the secure insertion of the **BENEFIT®** miniscrews. Implantation and appliance placement in one session.

Visit www.uniontech.it/easydriver for more information.

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
TADMAN®

The complete digital workflow realized with no efforts for the orthodontist. Based on the **BENEFIT®** System.

Visit www.tadman.de for more information.

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Specifications and design are subject to change without prior notice. All pictures, photos and product descriptions are the property of **psm** Medical GmbH. Utilization and copies by third parties have to be authorized by **psm** Medical GmbH. All rights reserved! All trademarks and brands mentioned in this brochure are subject to unrestricted regulations of the trademark laws valid and right of ownership of the registered owner.

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
Same instruments can be used for several **psm** products.
In order to find the right instruments for your use, they are indicated in colour.

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smartlock®

Introduction | Einführung

In collaboration with leading Orthodontic specialists, **psm** has developed the unique **mentoplate®** System for skeletal anchorage to facilitate maxillary advancement in early Class III treatment. The MENTOPLATE System offers a large variety of high precision titanium plates and screws for bone anchorage.

Especially in early Class III treatments in combination with RPE (**BENEFIT®** Hybrid Hyrax) the **mentoplates®** offer the possibility to treat children from the age of 8 years. Due to the fact that the fixation is performed in the sympheseal area, there is no need to wait for the eruption of the lower canines.

The plate insertion is made by means of a safe and stable fixation with selfdrilling miniscrews directly in the bone following a small incision. The maxillary protraction starts after finishing the RPE with intermaxillary elastic anchorage. The non-compliance treatment method without the need of a facemask shows excellent clinical results and is based on significant clinical research.

With the introduction of the **mentoplate®** System, **psm** has considerably enhanced the therapeutic diversity of early Class III treatment.

The new ANKA Plates used for maxillary or mandibular anchorage for a large variety of tooth movement applications with the added advantage that many **BENEFIT®** appliances can be attached to the plate giving amazing versatility.

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Literature notes | Literaturhinweise:

- Wilmes B., Nienkemper M., Ludwig B., et. al: Early Class III Treatment with a Hybrid Hyrax-Mentoplate Comination
J. Clin Orthod 45:1-7, 2011
- Wilmes B., Drescher D.: A miniscrew system with interchangeable abutments.
J. Clin Orthod 42:574-80; quiz 595, 2008
- Wilmes B., Drescher D., Nienkemper M.: A miniplate system for improved stability of skeletal anchorage.
J. Clin Orthod 43:494-501, 2009
- Ludwig B., Glasl B., Bowman J., Drescher D., Wilmes B. Miniscrew supported Class III Treatment with the Hybrid RPE Advancer.
J. Clin Orthod 44:533-539, 2010

The *mentoplate®* System | Das *mentoplate®* System

The Mentoplate-System provides a perfect solution for enhanced stability with special orthodontic bone anchored plates for intraoral fixation in the maxilla and the mandible. The System includes Titanium plates which are fixed to the bone surface with selfdrilling screws.

- Intrusion and extrusion of teeth
- Distal / mesial movement
- Space closure
- Maxillary protraction
- Maximum retention force
- Easy fixation with selfdrilling screws
- Reduced risk of tooth damage
- Precise control of tooth movement
- Ability to eliminate facemask

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





Abutments for Anka plate | Abutments für Anka Platte

33-54430 1ea. / 1St.

BENEFIT® Abutment Standard
BENEFIT® Abutment Standard

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

BENEFit® Schrauben im Spender 10 Stk.

Hybrid Hyrax® RPE for *BENEFit®* | Hybrid Hyrax® GNE for *BENEFit®*


The *Hyrax®* Ring and *Hyrax®* Abutment facilitate the construction and fixation of a *Hybrid Hyrax®* appliance. It allows welding or soldering with the *Hyrax®* Screw as well as little rotation movements when connecting the appliance to the *BENEFit®* Screws. The *Hyrax®* Ring is supplied with 1 fixation screw each.

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33-54463

1ea. / 1St.

BENEFit® Hyrax®, Abutment
BENEFit® Hyrax®, Abutment

bmx – *BENEFit®* Maxillary Xpander | bmx – *BENEFit®* Maxillary Xpander

After using the *BENEFit®* Screw for the *Hybrid Hyrax®* appliance very successfully in the last years the *bmx* now follows the Orthodontists request for an entirely *BENEFit®* mini-implant anchored appliance for Rapid maxillary expansion on young patients. Using the De Franco Screwguide (Page 18) the distance and the parallelity of the *BENEFit®* Screws can be exactly determined, directly on the patient. An impression is no longer necessary. Alternatively visualisation can of course be made by means of digital imaging or traditional impression taking.

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bmx, Gaumennahtspreizer

max. Expansion |
max. Ausdehnung:
10mm

Hybrid Hyrax® DIRECT | Hybrid Hyrax® DIRECT

for **BENEFIT®** DIRECT


33-64410

2ea. / 2 St.

Hyrax® Ring DIRECT 1-hole, stainless steel

Hyrax® Ring DIRECT 1-Loch, Implantatstahl

33-64414

1ea. / 1 St.

Hyrax® Ring DIRECT 1-hole with wire 1.1mm, 6cm, stainless steel

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33-64415

1ea. / 1 St.

Hyrax® Ring DIRECT 1-hole, median, with wire 1.1mm, 12cm, stainless steel

Hyrax® Ring DIRECT 1-Loch, median, mit Draht 1,1mm, 12cm, Implantatstahl

bmX DIRECT – **BENEFIT®** Maxillary Xpander

10

33-67406 6mm

33-67408 8mm

33-67410 10mm

bmX DIRECT, palatal expander

bmX DIRECT, Gaumennahtspreizer

max. Expansion |
max. Ausdehnung:

Plate distance
Plattenabstand

Video

12

33-67426 6mm

33-67428 8mm

33-67430 10mm

bmX DIRECT, palatal expander

bmX DIRECT, Gaumennahtspreizer

max. Expansion |
max. Ausdehnung:

For more details see the **psm** Smartlock catalogue

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
MARPE – Miniscrew Assisted Rapid Palatal Expansion

smartlock®

Video

Quadexpander

BENEFit® DIRECT Screw | BENEFit® DIRECT Schraube

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
33-64410

2ea. / 2 St.

Hyrax® Ring DIRECT 1-hole, stainless steel

Hyrax® Ring DIRECT 1-Loch, Implantatstahl


2X

ST

ST-33-64209-2	2.0 x 9mm
ST-33-64211-2	2.0 x 11mm

BENEFit® DIRECT Screws in a
package of 2 ea.

BENEFit® DIRECT Schrauben im
VPE 2 St.

33-64414

1ea. / 1 St.

Hyrax® Ring DIRECT 1-hole with wire 1.1mm, 6cm,
stainless steel

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10 Stück/Pkg.

Sterile packaging | Sterile Verpackungen

Implants – in single sterile packaging

The outer and inner packaging contains labels indicating UDI codelike article number, LOT number, expiration date, etc. for the patient folder and for full traceability. The sterile screws are contained in a small colour-coded screw carrier which protects the screw and at the same time acts as a pick up tray. This carrier is placed in a clear, sealed blister package. The screw can easily be removed from the carrier by means of the self retaining screwdriver.

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Universal Parallel guide 5-9mm
Universal Parallelführung 5-9mm

Video

Developed with
Entwickelt mit:
Dr. David J. De Franco
Milano, Italy

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The screwholder (33-10903) and drill (33-10914) can also be used with digital motion (page 34) and individually made insertion guides.

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33-825-22 45mm

TX 2.0 Screwholder for 10-61000, 1ea.

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BENEFit® DIRECT Screw | **BENEFit®** DIRECT Schraube

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
BENEFit® DIRECT TX Screw holder, 28 mm

BENEFit® DIRECT TX Schraubenhalter, 28 mm

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BENEFit® DIRECT Screw | **BENEFit®** DIRECT Schraube
TX Screw 2.0 | TX Schraube 2.0

33-18232 23 mm

33-18233 28 mm

33-18237 33 mm

Screwholder for contra-angled handpiece, 1ea.

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33-10903

1ea. / 1St.

BENEFit® / Digital Motion Screw holder with stop

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33-54515 30 mm

33-54518 50 mm

Screwdriver insert e.g. for abutment retaining screw,
BENEplate® fixation screws and Mobilizer activation,
1ea.

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Labor Schraubendreher, fest

1ea. / 1Stk.

Activating Screwdriver for mobilizer, not autoclavable,
acc. Dr. Michael Schubert, Regensburg, Germany

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INTRODUCTION | EINFÜHRUNG

Orthodontic Mini-Anchor System

The use of Mini-Anchor Screws or Temporary Anchorage Devices (TAD's) has become a standard of care within the orthodontic community over the past few years.

First offered in 2000 the **quattro®** TAD design has continuously been changed and improved to become the **quattro®** TAD family of today. The **quattro®** TAD's even after more than a decade of production are still unmatched in terms of versatility, quality and maximal success rate. With the **quattro®** MINI screw a new even smaller bracket head design is launched in 2012 offering a new dimension in patient comfort based on the outstanding **quattro®** design. Since year 2000 close to 100.000 **quattro®** TAD's have been placed in more than 70 countries worldwide.

The **quattro®** System is ideally enhanced with the unique **BENEFit/BENEplate** System (developed in cooperation with University of Düsseldorf, Germany) offering exchangeable abutments and connecting plates. This opens up new possibilities in palatal based appliances popularly used for molar distalization (see QR Video, page 26).

quattro® and **BENEFit®** a complete TAD System using the same basic instrumentation keeping inventory manageable and very cost effective.

Features:

- Clinically Proven Applications
- Titanium Implant Material
- Self-Drilling & Self-Tapping Technique
- Reliable & Secure Anchorage
- Innovative Head Attachment Design
- Stronger Conical Body Design
- Deeper Threads & Denser Pitch Design
- Various Screwhead Design
- Available in multiple Sizes

Benefits:

- Provides improved higher strength
- Immediate load-bearing of heavy forces
- Easy to combine with applications of orthodontic accessories
- Better mechanical retention
- Avoids gingival impingement
- Reduced treatment time
- Easy to insert and remove

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

quattro® PLUS

quattro® PLUS V


One Screw – Unlimited Possibilities

The *quattro®* TAD system offers 5 different screwhead designs which are perfectly adapted to the Orthodontist's appliances used daily. All *quattro* screws (1.3 mm to 2.3 mm) feature the identical and well-proven thread design for easy insertion and most secure and satisfactory treatment results.

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Based on the successful *quattro®* Screw the *quattro® MINI* Screw has been reduced in size in the head and neck area to meet compatible up with different anatomical conditions. The design of the selfdrilling thread is identical to the **psm quattro®** and **BENEFIT®** Screws which have demonstrated outstanding primary stability and very low failure rates for many years. The extremely sharp threads and the optimized thread pitch facilitate the insertion tremendously and make the *quattro® MINI* an ideal addition to your selection of Standard Mini Anchor Screws.

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

ST

All implants are supplied in a pre-sterilized single package

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RH: Round head design | Runder Kopf Design


BH: Bracket head design | Bracket Kopf Design

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


Emergency Screw, 1 ea.
Notfallschraube, 1 Stk.


quattro® HOOK

Ø1.5

Size (Ø x length)
Abmessung (Ø x Länge)

ST

ST-33-21507	1.5 x 7 mm
ST-33-21509	1.5 x 9 mm
ST-33-21511	1.5 x 11 mm


Ø2.0

Size (Ø x length)
Abmessung (Ø x Länge)

ST


ST-33-22007	2.0 x 7 mm
ST-33-22009	2.0 x 9 mm
ST-33-22011	2.0 x 11 mm

quattro® HOOK Screw, Titanium, 1 ea.
quattro® HOOK Schraube, Titan, 1 Stk.

ST-33-22109	2.3 x 9 mm
--------------------	------------

Emergency Screw, 1 ea.
Notfallschraube, 1 Stk.

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





DENTAL attachment
DENTAL Ansatz

12-12625

Drill, 1.0x33mm WL 15mm, green, for 1.8mm screws
Bohrer, 1,0x33mm AT 15mm, grün, für 1,8mm Schrauben

10-12208

Drill, 1.2x 22mm WL 8.5mm, green, for 1.8mm screws
Bohrer, 1,2x 22mm AT 8,5mm, grün, für 1,8mm Schrauben

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The extension of 16 mm works for all rotary instruments of the **quattro®** and **BENEfit®**-System with a latch type DENTAL attachment (max. rpm 600).

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
Drill stop | Bohrerstop

The drill stop prevents the clinician from unintentionally drilling too deep. They assure an accurate working length from 7 to 13 mm in 2 mm increments. By sliding the stop sleeve over the appropriate **psm** drill, the working length is adjusted to the desired depth. Removal of the sleeve is easy by using the buffer.

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All *quattro®* screwholders are marked with an arrow indicating the slot, hole or hook position on the *quattro®* screw during insertion.

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
Cordless Battery Screwdriver with Torque Calibration System for the fixation of **quattro®** and **BENEFIT®** Miniscrews.

The iSD 900 will guarantee precise torque control due to the unique – easy to use – torque calibration system (10 – 40 Ncm). Time-saving and precise, riskfree insertion and tightening of screws and abutments. The complete Set contains a reliable and user friendly handpiece for Dental attachments. Reverse mode is indicated with sound. Standard rechargeable AAA Ni-Mh batteries can be used.

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**33-18240**

1ea. / 1Stk.

Torque ratchet, handle with insert 33-18234 for
quattro® / BENEFit®

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■ ■ ■
73-31960

1ea. / 1Stk.

quattro® large tray for instruments and implants, without content, with lid

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33-11500		quattro® / BENEFit® Instrument Set consisting of:	quattro® / BENEFit® Instrumenten Set bestehend aus:
		Storage	Aufbewahrung
73-31960	Instruments tray	Instrumenten-Behälter	
		Instruments	Instrumente
10-61000	Screwdriver handle	Schraubendreher Handgriff	
33-18230	Screwholder QB, long, for 10-61000	Schraubenhalter QB, lang, für 10-61000	
33-825-22	BENEFit® DIRECT TX Screwholder for 10-61000	BENEFit® DIRECT TX Schraubenh. f. 10-61000	
33-18236	Thumb screw mandrel for Dental shaft instr.	Handgriff für Instrumente mit Dentalansatz	
33-18232	Screwholder QB, short, 23 mm, DENTAL	Schraubenhalter QB, kurz, 23 mm, DENTAL	
33-825-33	BENEFit® DIRECT TX Screwholder, short	BENEFit® DIRECT TX Schraubenhalter, kurz	
33-54512	BENEFit® Screwholder for Abutment, 20mm	BENEFit® Eindrehinstr. für Abutment, 20mm	
12-12625	Drill, 1.0 x 33 mm, WL 15 mm, DENTAL, green	Bohrer, 1,0 x 33 mm AT 15 mm, DENTAL, grün	
10-67513	Drill, 1.4 x 33 mm, WL 15 mm, DENTAL, red	Bohrer, 1,4 x 33 mm AT 15 mm, DENTAL, rot	
33-18288	Laboratory Screwdriver for Abutments etc.	Labor Schraubendreher für Abutments etc.	
33-18244	Drill extension +16 mm, DENTAL	Bohrer Verlängerung +16 mm, DENTAL	

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73-31971

1ea. / 1Stk.

quattro® small tray for instruments, without content, with lid

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Ø 2.0 mm screw

Catalogues | Kataloge

Even more innovative **psm** products are presented in our other catalogues. Please request it from your local **psm** distribution partner or download from our website www.psm-medical.com.

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BENEFit

With the introduction of the **BENEFit®** System **psm** has considerably enhanced the therapeutic diversity of the treatments available for Mesialization, Distalization, RPE as well as uprighting or anchorage of the dentition. With **BENEFit®**, **psm** was the first company to offer a full range of skeletal anchorage appliances for the palate – a trend that has become „state of the art“ in Orthodontics.

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smartlock

The new **psm – smartlock®** system is largely based on the successful and proven **BENEFit®** system. **smartlock®** focuses on the essentials: Unique comfort as well as safe and easy handling that doctors and patients can rely on. The trend-setting concept, combined with the best materials and the most modern manufacturing processes, guarantees a treatment result at the highest level.

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mentoplate

The Mentoplate-System provides a perfect solution for enhanced stability with special orthodontic bone anchored plates for class III correction. The System includes Titanium plates which are fixed to the bony surface with selfdrilling screws.

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Digital skeletal anchorage | Digitale skelettale Verankerung

A new and innovative method for the secure insertion of the **BENEFit®** miniscrews.
Implantation and appliance placement in one session.

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BENEFit

For more information please refer to our **BENEFit®** System – Handout or our website www.psm-medical.com

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
Introduction | Einführung

Difference **BENEFit®** & **BENEFit® DIRECT** | Unterschied **BENEFit®** & **BENEFit® DIRECT**

BENEFit

IMPLANT FIRST

The traditional workflow to insert the classic **BENEFit®** Implants first and then build the appliance on top has its proven advantages when it comes to distalization (Beneslider) and mesialization (Mesialslider) as there are no points to fix the appliance in the palate before screw insertion.

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BENEFit DIRECT

BENEFit® DIRECT

The new **BENEFit® DIRECT** Screws show their big advantages when an APPLIANCE FIRST workflow is preferred. Without any additional software the appliance (i.e. Hybrid Hyrax, Marpe) can be prepared on the patient model with our laboratory tools and pre-welded **bmX** DIRECT Expanders. The appliance can be fixed on the molars and the **BENEFit® DIRECT** Screws secure the appliance safely in the desired position. The threaded head design ensures a secure and firm connection between the DIRECT plates and the implant without any tilting.

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Products Produkte

BENEFit® DIRECT Screws | BENEFit® DIRECT Schrauben


Ø2.0

Size (Ø x length)
Abmessung (Ø x Länge)

ST■

ST-33-64207	2.0 x 7mm
ST-33-64209	2.0 x 9mm
ST-33-64211	2.0 x 11mm
ST-33-64213	2.0 x 13mm

BENEFit® DIRECT Screw 2.0 mm

BENEFit® DIRECT Schraube 2,0 mm

2X

ST■

ST-33-64209-2 2.0 x 9mm

ST-33-64211-2 2.0 x 11mm

BENEFit® DIRECT Screws in a package of 2 ea.

BENEFit® DIRECT Schrauben im VPE 2 St.

NOTE | HINWEIS:




ST■ All implants are supplied in sterile single use packaging.

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median

paramedian

Hybrid Hyrax DIRECT | Hybrid Hyrax DIRECT

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

1ea. / 1 St.

33-64415

Hyrax® Ring DIRECT 1-hole, median, with wire 1.1mm, 12 m, stainless steel

Hyrax® Ring DIRECT 1-Loch, median, mit Draht 1,1mm, 12cm, Implantatstahl

1ea. / 1 St.

bmx DIRECT – BENEFit® Maxillary Xpander

10

33-67406 6mm

max. Expansion |

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33-67408 8mm

33-67410 10mm

bmx DIRECT, palatal expander

bmx DIRECT, Gaumennahtspreizer

Plate distance

Plattenabstand

Video

12

33-67426 6mm

max. Expansion |

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33-67428 8mm

33-67430 10mm

bmx DIRECT, palatal expander

bmx DIRECT, Gaumennahtspreizer

Micro glue sticks can be found on page 20.

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2ea. / 2 St.


**33-64879**

smartlock® Pin Set
smartlock® Pin Set

2ea. / 2 St.

**33-65000**

Laboratory Instrument Set
Labor Instrumenten Set

 **Video**

Laboratory Products | Labor Produkte

33-65000	Laboratory Instrument Set consisting of:	Labor Instrumenten Set bestehend aus:
Instruments		Instrumente
33-18236	Thumb screw for Dental mandrel	Handgriff für Klingen mit Dentalansatz
33-825-33	BENEFit® DIRECT TX Screwholder	BENEFit® DIRECT TX Schraubenthalter
33-54512	Screwholder for fixation screw, 20 mm	Schraubenthalter für Fixierschraube, 20 mm
33-18232	Screwholder for Smartlock Cone Pins 33-64880	Schraubenthalter für Smartlock Kegelpin 33-64880

Instruments Instrumente

Caption | Legende

The instruments can be used for several **psm** products. In order to find the right instruments for your use, they are indicated in colour.

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The extension of 16 mm works for all rotary instruments of the quattro® and **BENEFit®**-System with a latch type DENTAL attachment (max. rpm 600).

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33-54515 30 mm

33-54518 50 mm

Screwdriver insert e.g. for laboratory and fixation screws, 1ea.

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33-18288

1ea. / 1St.

Laboratory Finger Screwdriver, fixed

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Cordless Battery Screwdriver with Torque Calibration System for the contolled fixation of **psm** Miniscrews.

The iSD 900 will guarantee precise torque control due to the unique – easy to use – torque calibration system (10 – 40 Ncm). Time-saving and precise, riskfree insertion and tightening of screws and abutments. The complete Set contains a reliable and user friendly handpiece for Dental attachments. Reverse mode is indicated with sound. Standard rechargeable AAA Ni-Mh batteries can be used.

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

33-64890

64 ea. / 64 St.

Micro Stix, glue sticks
Micro Stix, Klebesticks

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* as long as stock lasts
* solange Vorrat reicht

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The outer and inner packaging contains labels indicating article number, LOT number, expiration date, etc. for the patient folder and for full traceability.

The sterile screws are contained in a small color-coded screw carrier which protects the screw and at the same time acts as a pick up tray. This carrier is placed in a clear, sealed blister package. The screw can easily be removed from the carrier by means of the self retaining screwdriver.

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Ø 2.0 mm screw

Catalogue | Kataloge

Even more innovative **psm** products are presented in our other catalogues. Please request it from your local **psm** distribution partner or download from our website www.psm.ms.

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quattro

The **quattro®** Mini Anchor System is known as one of the first and most successful orthodontic Mini-Implant Systems worldwide. Based on this system the **BENEFit®** screws were developed. This enables the user to utilize the entire scope of orthodontic implantation with unmatched possibilities and compatible instrumentation.

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EASY DRIVER[®] THE GUIDING SYSTEM FOR TADS

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BENEFit

For more information please refer to our **BENEFit®** System – Handout or our website www.psm.ms

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PSM North America, Inc.
PSM Medical GmbH